

Quantitative Analysis of Meat Spoilage using VIS/NIR Spectral Imaging

J.M. Carstensen,^{1*} E. Panagou² and G-J.E. Nychas²

¹Videometer A/S, Lyngso Allé 3, 2970 Horsholm, Denmark, E-mail: jmc@videometer.com

²Agricultural University of Athens, Department Food Science, Technology & Human Nutrition, Lab of Microbiology & Biotechnology of Foods, Iera Odos 75, Athens 11855 Greece

Carried out with the financial support of the Commission of the European Communities, specifically the RTD programme, IP Project No. 211638 SYMBIOSIS-EU

Meat spoilage study



- Minced beef shelf life
- Stored at 5° under
 - Aerobic conditions (normal atmosphere)
 - Modified atmosphere packaging (MAP) 40% CO₂; 30% O₂; 30% N₂
- Measured at 12 hour intervals over 7 days

Data acquisition and microbiological analysis: Ammor, M.S., Argyri, A., Nychas, G.-J.E. (2009) Rapid monitoring of the spoilage of minced beef stored under conventionally and active packaging conditions using Fourier transform infrared spectroscopy in tandem with chemometrics. *Meat Science* 81, 507-514



Quantification of meat spoilage


- 50+ methods have been reported and are reasonably documented
- There is still a need for a method that is
 - Fast
 - Non-contact (repetitive measurements on same sample)
 - Non-labor-intensive
 - Cost-efficient
 - Deals effectively with the heterogeneity of the sample
- Spectral imaging is a serious candidate



Measurements performed

- Microbiological analysis
 - Total viable count (TVC)
 - Pseudomonas
 - Brochothrix thermosphacta
 - Lactic acid bacteria (LAB)
 - Enterobacteriaceae
 - Yeast and moulds
- Other analyses
 - pH, sensory quality
- Multispectral imaging (Videometer)






Laboratory device for multispectral imaging



- Up to 20 spectral bands in the range 360 nm to 1050 nm
- Up to 1600×1200 pixels per band
- Very homogeneous and diffuse illumination
- Strobed LED light source



Meat spoilage effects in VIS-NIR

- Change existing muscle pigment states
- Generate new pigments/metabolites
- Change the surface chemical abundance on relevant compounds



Muscle pigments: myoglobins




FIG. 3. Absorption curves for Mb, MbCO, and MMbCN in the visible region

- Myoglobin / Oxymyoglobin
 - It is the primary oxygen-carrying /storing pigment of muscle tissues
- Metmyoglobin
 - is the oxidized form of myoglobin.

Source: W. Bowen, 1948.



OxyMb



MetMb



DeoxyMb

Image source Danish Meat Research



Before storage


1a



1b



1c




1d



Aerobic storage

A4a




A8a



A12a




A13a




MAP storage

M4a




M8a







M12a






M13a



Microbiological analysis





Canonical discriminants

Find the \mathbf{a} that maximizes the Rayleigh quotient


$$R(\mathbf{a}) = \frac{\mathbf{a}^T \boldsymbol{\Sigma}_S \mathbf{a}}{\mathbf{a}^T \boldsymbol{\Sigma}_N \mathbf{a}}$$

- where $\boldsymbol{\Sigma}_S = \mathbf{A}$, the among class dispersion, and $\boldsymbol{\Sigma}_N = \mathbf{W}$, the within class dispersion
- Maximizing this quotient originates back to Fischer (1936)
- Eigenvectors are called canonical discriminants functions or simply canonical discriminants
- The number of canonical discriminants CDs is limited by the rank of \mathbf{A} . With two classes there is only one CD. With three classes there are two CDs
- Typically a supervised technique




Canonical discriminant images


1a




1b



M13a




A13a



Canonical discriminant histograms

Histograms




- Red – Before storage (1a)
- Brown – Before storage (1b)
- Cyan – MAP storage (M13a)
- Blue – Aerobic storage (A13a)




Canonical discriminant plots

90% quantile CDF (negated)




MNF 1






Upper left: No storage

Upper right: Aerobic

Lower right: MAP








Upper left: No storage

Upper right: Aerobic

Lower right: MAP







Upper left: No storage

Upper right: Aerobic

Lower right: MAP






Upper left: No storage

Upper right: Aerobic

Lower right: MAP



Conclusions

- Color and surface chemistry changes during meat spoilage may be monitored using VideometerLab
- Heterogeneity of changes becomes apparent and may be measured
- A canonical discriminant function (CDF) shows that
 - Aerobic storage gives a large a bimodal change in color
 - Aerobic storage gives most color change from 50 to 90 hours storage
 - MAP storage gives a smaller but significant and unimodal change in color
 - The changes above can be directly quantified from the histogram of the CDF

